


Using NI Elvis and MyDAQ to Support Remote Labs for STEM Education

8/1/24 AGBell



Using NI Elvis and MyDAQ to Support Remote Labs for STEM Education

This presentation will discuss the detailed designs using NI MyDAQs and NI ELVIS IIIs to build circuits that can be remotely accessed from outside of the college. Lab experiments with RC and RL circuits, relays, MyDAQs, and the NI ELVISmx software suite can be used to evaluate RC and RL circuits. Likewise, lab experiments for active filters can be evaluated using NI ELVIS IIIs and NI Measurement Live software. Both the MyDAQ and ELVIS III boxes connect to computers via a USB interface. Each of the four rack-mounted computers is connected to both an ELVIS III and MyDAQ.






How to use GRE computers and associated equipment using Beyond Trust

- A. To access beyond go to <https://access.ivytech.edu/login/login>
- B. Select Use SAML Authentication hyperlink under the Login button
- C. Enter your Ivy Tech credentials
- D. Select Launch Privileged Web Access Console
- E. Under the Jump Items column highlight FortWayne-GRE-1337 (4 servers should appear numbering 002-005)
- F. Select which server you would like to use. (Depending on slots available)
Select the orange JUMP button

We use a common username and password for the 4 GRE computers.

Due the cost of Beyond Trust we will be using a special VPN on the Ivy Tech this fall using Cisco AnyConnect. Students will still a common username and password, but they will need to load the Cisco AnyConnect software on their own computers ...



How to use GRE computers and associated equipment using Beyond Trust (cont)

G. The windows login screen should appear.

H. Treat this window like Windows desktop pick the student account and enter the common password

II. Setup MyDAQ

A. National Instruments software suite

1. Navigate to windows icon in the left bottom corner
2. Find National Instruments drop down menu
3. Locate file shortcut NI ELVISmx Instrument Launcher and right click to start




B. Using the NI ELVISmx Instrument Launcher



II. Setup MyDAQ (cont)

B. Using the NI ELVISmx Instrument Launcher

1. For an oscilloscope experiment start the following instruments: Oscilloscope >> Function generator >> Digital Writer
2. Schematic for the relay logic used for EECT series of experiments. AC input signal is diverted using rly0 to either feed capacitor or inductor side of circuit. The relays 1-6 are used to select the sample being tested.




II. Setup MyDAQ (cont)

B. Using the NI ELVISmx Instrument Launcher

3. Function Generator operation type in frequency, Amplitude, then the Run button.
(The frequency should appear)

4. Oscilloscope first enable channels A10 and A11 (with nothing selected on digital writer nothing selected both channels should look very similar in frequency and Vpp)
Then press the Run button

5. Digital Writer open the Lines to Write drop down menu select 0 – 7



III. Conduct Experiment

B. Using the NI ELVISmx Instrument Launcher

- Lay out is bit0 switches the function generator between the capacitors and the inductors. Bit1-3 selects capacitors, and bit4-6 selects the inductors to be tested
- Each selector bit corresponds to the relay (bit1 = relay1, bit2 = relay2 ... etc.) After selecting the desired bits press the Run button. You may switch between them
- With the selector bit1 activated a 0.47 μ f and a 1k ohm resistor is being displayed (the signal on AI1(Blue) reflects the change). Take the measurement from channel 1 and compare to the results of the calculations. Note: The nominal value written on component may differ from the actual value.



III. Conduct Experiment (cont)

B. Using the NI ELVISmx Instrument Launcher

d. DON'T FORGET TO DEACTIVATE THE BIT BEFORE MOVING ON TO THE NEXT EXPERIMENT!!!

e. After Deactivating the bit, move on to the next experiment

f. To access the inductor side of the experiment first select bit0 (this engages the switching relay), then select an inductor bit4-6

g. Performing the inductor set of experiments. Bit0 and bit4-6 is actuated, with 1 kHz applied at 0.5Vpp.



III. Conduct Experiment (cont)

B. Using the NI ELVISmx Instrument Launcher

6. AC Frequency Response - Bode plots

- a. Within the NI ELVISmx Instrument Launcher window Bode Analyzer and Digital Writer.
- b. Control the relay bits the same way as in previous experiments apply the logic needed to examine the needed components
- c. Set the frequency, turn the cursors on and run this may take a while to finish and the Run button
- d. Set the next experiment logic up and repeat

Stop all running windows and log out

Results need to be saved to your Google Drive!



Key feature

Pros

Very flexibility

Software and hardware are mature.

Cost is low and students can buy their own MyDAQ.

Could support active and passive circuits.

USB interface

Cons

Bandwidth is not real high (20 KHz)



No LCR function


May be a little trick for some students to use.



NI MyDAQ

Ivy Tech designed, and 3D printed a 7-relay holder using Solidworks that supported the used of a full-sized solderless breadboard and seven 5V relays.





Frequently Used Jump Items

REFRESH

CFTI-1337-005

Jump Client

FortWayne-GRE-1337

Active [ON]

Last accessed at: 03/16/2023 1:43 AM

CFTI-1337-004

Jump Client

FortWayne-GRE-1337

Active [ON]

Last accessed at: 04/27/2023 9:40 AM

CFTI-1337-003

Jump Client

FortWayne-GRE-1337

Active [ON]

Last accessed at: 04/27/2023 9:24 AM

CFTB-1337-N01

Jump Client

FortWayne-CFTB-1337


Active [ON]

Last accessed at: 04/25/2023 6:32 PM


FortWayne-GRE-1337

Name ▲	Jump Method	Group	Status	Last Accessed		
+	CFTI-1337-002	Jump Client	FortWayne-GRE-1337	Active [ON]	04/27/2023 9:12 AM	JUMP
+	CFTI-1337-003	Jump Client	FortWayne-GRE-1337	Active [ON]	04/27/2023 9:24 AM	JUMP
+	CFTI-1337-004	Jump Client	FortWayne-GRE-1337	Active [ON]	04/27/2023 9:40 AM	JUMP
+	CFTI-1337-005	Jump Client	FortWayne-GRE-1337	Active [ON]	03/16/2023 1:43 AM	JUMP

Launch BeyondTrust and select one of the FortWayne-GRE-1337 computers. If the computers are not online contact Professor Bell



Open Measurements Live



Measurements Live
https://measurementslive.ni.com
Apps IvyTech.edu MyIvy IvyTechEngineering Slimjet search Fastest VPN in the... Slimjet homepage
Resources

MeasurementsLive

Select Measure. You are more than welcome to review the other options.



FIRST TIME HERE?
Let's get you set up. Follow these steps to get ready to use the NI ELVIS III device.



MEASURE
Already set up? Connect your device and use the instruments.



DEVICE SIMULATION
Learn to use the instruments before connecting to a real device.

© 2019 National Instruments | Legal | Privacy | Export | Data Use Notice | Notices
NATIONAL INSTRUMENTS

Measurements Live
https://measurementslive.ni.com/measure.html
Apps IvyTech.edu MyIvy IvyTechEngineering Slimjet search Fastest VPN in the... Slimjet homepage
Resources

No device connected **Manage device connection**

USB Network

Connect the device's rear USB-C port to your computer. Click Connect to pair the device.



Warning: A popup window may appear.

Your name here

Need the software?

Connect

Under the USB option
Select the "Connect"

Click here
to add instruments

Measurements Live
https://measurementslive.ni.com/measure.html
NI-ELVISIII-0315ACA1 Manage device connection
NI-ELVISIII-0315ACA1
Device: NI ELVIS III
Application Board: ---
Instruments
Device Details
Disconnect
Resources
Click here to add instruments
Select the "add instruments" option

Measurements Live
https://measurementslive.ni.com/measure.html
NI-ELVISIII-0315ACA1 Manage device connection

Select the Variable Power Supply and Bode Analyzer

- Oscilloscope [?] Available
- Function and Arbitrary Waveform Generator [?] Available
- Data Generator [?] Available
- Digital Multimeter [?] Available
- Variable Power Supply [?] Available
- Bode Analyzer [?] Available
- Current-Voltage Analyzer [?] Available
- Logic Analyzer and Pattern Generator [?] Available Separate window
- Digital I/O [?]

10:35 AM 4/27/2023

Measurements Live
https://measurementslive.ni.com/measure.html
NI-ELVISIII-0315ACA1 Manage device connection
VPS Run
Positive (+) Static
Voltage - - - V
Current - - - A
Power - - - W
1 V, 500 mA
Voltage [?] 1V
Current limit [?] 500 mA
Negative (-) Inactive

Adjust the Positive voltage to 9V
Select the Negative Voltage dropdown
and pick "Statics". The set the
Negative voltage to -9V

Measurements Live
https://measurementslive.ni.com/measure.html
NI-ELVISIII-0315ACA1 Manage device connection
VPS Run
Power - - - W
9 V, 500 mA
Voltage [?] 9 V
Current limit [?] 500 mA
Negative (-) Static
Voltage - - - V
Current - - - A
Power - - - W
-9 V, -500 mA
Voltage [?] -9 V
Current limit [?] -500 mA

The voltages should look like what is shown in the left.

You can adjust the current limits but this is not absolutely required. +/- 10 mA would be fine or just leave it alone.

Measurements Live
https://measurementslive.ni.com/measure.html
NI-ELVISIII-0315ACA1 Manage device connection

Bode Analyzer Run

Frequency Cu: Off

Gain (dB)
Phase (°)

Power - - - W

9 V, 500 mA
Voltage [?] 9 V
Current limit [?] 500 mA

Negative (-) Static

Voltage - - - V
Current - - - A
Power - - - W

-9 V, -500 mA
Voltage [?] -9 V
Current limit [?] -500 mA

Type here to search 55°F 11:11 AM 4/27/2023

Select the Bode Analyzer also and it will launch.

Measurements Live
https://measurementslive.ni.com/measure.html
NI-ELVISIII-0315ACA1 Manage device connection

Bode Analyzer

Run

Frequency | Cu: Off

Stimulus channel

Start frequency	Stop frequency
10 Hz	1 MHz
Steps per decade	Peak amplitude
10	1 V

FGen/Arb CH1 - Oscilloscope CH1 [?]

Response channels

Response 1	<input checked="" type="checkbox"/>
Oscilloscope CH2	<input type="checkbox"/>
Response 2	<input type="checkbox"/>
Response 3	<input type="checkbox"/>

Reference channels

Reference 1	<input type="checkbox"/>
-------------	--------------------------

Voltage [?] -9 V

Current limit [?] -500 mA

Windows taskbar: Type here to search, 55°F, 11:12 AM 4/27/2023

If you select the little “gear” icon you can adjust the Start and Stop frequency (10 Hz to 1 MHz could be changed to 100 Hz to 100KHz), number of Steps per decade (if you want more resolution), Peak Amplitude (1V should be fine) and the response you want to observe.

Measurements Live
https://measurementslive.ni.com/measure.html
Apps IvyTech.edu MyIvy IvyTechEngineering Slimjet search Fastest VPN in the... Slimjet homepage

MeasurementsLive

NI-ELVISIII-0315ACA1 Manage device connection

Bode Analyzer

Run

Frequency | Cu: Off

Gain (dB)

Phase (°)

Stimulus channel

Start frequency: 10 Hz | Stop frequency: 1 MHz

Steps per decade: 10 | Peak amplitude: 1 V

FGen/Arb CH1 - Oscilloscope CH1

Response channels

- Response 1
- Oscilloscope CH2
- Response 2
- Response 3

Reference channels

- Reference 1

Current limit [?] 500 mA

Negative (-) Static

Voltage - - - V

Current - - - A

Power - - - W

-9 V, -500 mA

Voltage [?] -9 V

Current limit [?] -500 mA

Windows taskbar: Type here to search, 55°F, 11:12 AM 4/27/2023

Response 1 should be the Low-Pass filter response. Response 2 and 3 will be the High-Pass and Band-Pass responses.

MeasurementsLive
NI-ELVISIII-0315ACA1 Manage device connection

Bode Analyzer Run

Frequency | Cu: Off

Gain (dB)

Phase (°)

Stimulus channel

Start frequency	Stop frequency
10 Hz	1 MHz
Steps per decade	Peak amplitude
10	1 V

Response channels

Response 1	<input checked="" type="checkbox"/>
Oscilloscope CH2	<input type="checkbox"/>
Response 2	<input type="checkbox"/>
Response 3	<input type="checkbox"/>

Reference channels

Reference 1	<input type="checkbox"/>
-------------	--------------------------

VPS Run

Power

Voltage

Current

Power

-9 V, -500 mA

Voltage [-?] -9 V

Current limit [-?] -500 mA

To start the test click the "run" button on the Variable Power Supply *first* and then the "run" button on the Bode Analyzer *second*.

Measurements Live
https://measurementslive.ni.com/measure.html
NI-ELVISIII-0315ACA1 Manage device connection

Bode Analyzer

Run

Frequency | Cu: Off

Gain (dB)

Phase (°)

10 100 1k 10k 100k 1M

Negative (-) Static

Voltage	-9.0	V
Current	-3.3	mA
Power	29.5	mW

▼ -9 V, -500 mA

Type here to search | 55°F | 11:19 AM 4/27/2023

The filter response looks like a Lowpass filter. We can add cursors to measure the filter characteristics, take a picture and export the data as a CSV file for further analysis.

Measurements Live
https://measurementslive.ni.com/measure.html
Apps IvyTech.edu MyIvy IvyTechEngineering Slimjet search Fastest VPN in the... Slimjet homepage
MeasurementsLive Resources
NI-ELVISIII-0315ACA1 Manage device connection
Bode Analyzer Run

Frequency Cu: Manual

	C1 Response 1	C2 Response 1	Δ Value
Frequency	10 Hz	709.53 Hz	699.53 Hz
Gain	0.01229 dB	-3.0546 dB	3.0669 dB
Phase	-0.88144°	-45.415°	44.534°

Gain (dB)
Phase (°)

Negative (-) Static
Voltage -9.0 V
Current -3.5 mA
Power 31.1 mW
-9 V, -500 mA

To add the cursors select the dropdown arrow next to Cu: Off and change to Manual. You will also need to unlock them by selecting the little "lock" icon at the top of the Bode Analyzer window.

Measurements Live
https://measurementslive.ni.com/measure.html
NI-ELVISIII-0315ACA1 Manage device connection

Bode Analyzer

Run

Frequency | Cu: Manual

	C1 Response 1	C2 Response 1	\Delta Value
Frequency	10 Hz	709.53 Hz	699.53 Hz
Gain	0.01229 dB	-3.0546 dB	3.0669 dB
Phase	-0.88144°	-45.415°	44.534°

Gain (dB)

Phase (°)

9 V, 500 mA
Voltage [?] 9V
Current limit [?] 500mA
Negative (-) Static

Voltage -9.0 V
Current -3.5 mA
Power 31.1 mW

-9 V, -500 mA

Type here to search | 55°F | 11:30 AM 4/27/2023

This Low-Pass filter has a 3dB point at about 700 Hz.

The screenshot shows the MeasurementsLive web application interface. The main window displays a Bode Analyzer for a device labeled NI-ELVISIII-0315ACA1. The interface includes a table of response data, a Bode plot, and various configuration panels. A 'Take Screenshot' dialog box is open on the right side, with a red circle highlighting the camera icon in the top right corner of the application. A text box explains that selecting this icon allows for capturing the current response data, including cursor positions and values.

Frequency	C1 Response 1	C2 Response
Frequency	10 Hz	709.53 Hz
Gain	0.0108 dB	-3.0559 dB
Phase	-0.79406°	-45.274°

Stimulus channel

Start frequency: 10 Hz, Stop frequency: 1 MHz
Steps per decade: 10, Peak amplitude: 1 V

Response channels

- Response 1:
- Oscilloscope CH2:
- Response 2:
- Response 3:

Take Screenshot Dialog:

Choose instrument(s) to capture

- Bode Analyzer
- Variable Power Supply

Take Screenshot

By selecting the camera icon you can take a picture of the response to include the cursor locations and values.

Recycle bin
Microsoft Edge
Google Chrome
Measure Live
NI ELVIS Instrum
Briefca

Channel	Response 1	Response 2	Δ Value
Frequency	19.70	799.53 70	499.53 70
Gain	8.8188 dB	-3.8533 dB	2.8647 dB
Phase	-8.79496°	-45.274°	44.48°

Gain (dB)
Phase (°)

Stimulus Channel

Start Frequency	10 Hz
Stop Frequency	1 MHz
Steps Per Decade	10
Peak Amplitude	1 V

Response Channels

Response Name	State	Channel Label
Response 1	Enabled	Oscilloscope CH2
Response 2	Disabled	
Response 3	Disabled	

Reference Channels

Reference Name	State	Mode	Source File Name	Channel
Reference 1	Disabled			
Reference 2	Disabled			

4/27/2023 14:37:09

Picture of the response can be saved to your Google Drive for later use.

MeasurementsLive
NI-ELVISIII-0315ACA1 Manage device connection

Bode Analyzer

Run

Frequency	C1 Response 1	C2 Respo
Frequency	10 Hz	709.53 H
Gain	0.0108 dB	-3.0559
Phase	-0.79406°	-45.274°

Stimulus channel
Start frequency: 10 Hz, Stop frequency: 1 MHz
Steps per decade: 10, Peak amplitude: 1 V

Response channels
Response 1:
Oscilloscope CH2:
Response 2:
Response 3:

Choose data to export
Bode Analyzer
 CSV
Export Data

bode_20230427_1....png

2:42 PM 4/27/2023

Likewise the data can be exported to a CSV file for more detailed analysis and saved in your Google Drive.

Filter Response - Google Drive x +

https://drive.google.com/drive/folders/196Ky_s5RpQHwE8DMvhwfNhnjsjkYCMtN

Apps IvyTech.edu MyIvy IvyTechEngineering Slimjet search Fastest VPN in the... Slimjet homepage

Drive

Search in Drive


+ New

- Priority
- My Drive
- Shared drives
- Shared with me
- Recent
- Starred
- Trash
- Storage
400.35 GB used

M... > Filter R... ▾

Files Name ↑

bode_20230427_143740.png



bode_20230427_144427.csv

Who has access

- A

2 uploads complete

- bode_20230427_144427.csv ✓
- bode_20230427_143740.png ✓

Type here to search

62°F 2:49 PM 4/27/2023

Save your work in you Google Drive via MyIvy

Measurements Live
https://measurementslive.ni.com/measure.html
Apps IvyTech.edu MyIvy IvyTechEngineering Slimjet search Fastest VPN in the... Slimjet homepage

MeasurementsLive

NI-ELVISIII-0315ACA1 Manage device connection

Bode Analyzer

Run

Frequency	C1 Response 1	C2 Respo
Frequency	10 Hz	709.53 H
Gain	0.0108 dB	-3.0559
Phase	-0.79406°	-45.274°

Stimulus channel
Response channels
Response 1
Oscilloscope CH2
Response 2
Response 3
Reference channels
Reference 1
Reference 2

9 V, 500 mA
Voltage [?] 9 V
Current limit [?] 500 mA
Negative (-) Static
Voltage -9.0 V
Current -3.7 mA

bode_20230427_14...xls bode_20230427_14...xls Show all

Type here to search S&P... 2:53 PM 4/27/2023

Select Response 2 and then Response 3 to get the High-Pass and Band-Pass responses, respectively.

Measurements Live

https://measurementslive.ni.com/measure.html

NI-ELVISIII-0315ACA1 Manage device connection

Bode Analyzer

Run

Frequency	1.0061 kHz	10.351 kHz	9.3451 kHz
Gain	-3.0107 dB	-0.01201 dB	2.9987 dB
Phase	-269.71°	-352.69°	82.972°

Gain (dB)

Phase (°)

C1 C2

Voltage 9.0 V
Current 2.9 mA
Power 25.7 mW

9 V, 500 mA

Voltage [?] 9V

Current limit [?] 500mA

Negative (-) Static

Voltage -9.0 V
Current -3.5 mA

bode_20230427_14....xls

3:14 PM 4/27/2023

High-Pass response with a 3dB point @ 1.0KHz

Measurements Live
https://measurementslive.ni.com/measure.html
NI-ELVISIII-0315ACA1 Manage device connection

Bode Analyzer

Frequency	985.89 Hz	2.3037 kHz	1.3179 kHz
Gain	4.6725 dB	-19.219 dB	23.892 dB
Phase	-170.7°	-266.9°	96.196°

Gain (dB) vs Frequency (Hz) and Phase (°) vs Frequency (Hz) plot. The plot shows a resonance peak at 985.89 Hz. Vertical dashed lines mark C1 and C2.

Run

Band-Pass center frequency @ 985Hz

9.0 V
2.7 mA
24.7 mW

9 V
500 mA

Static

-9.0 V
-3.5 mA

bode_20230427_14...xls

Questions



Andy Bell
Department Chair – Engineering
Ivy Tech Community College – Northeast
Phone: 260-481-2288 : Fax: 260-480-2052
abell118@ivytech.edu
SDKB Technology Center, Room TC1240R,
3800 N. Anthony Blvd.,
Fort Wayne, IN 46805
SL avatar = ivytechengineer

